



SRT-LVL

LVV SIZE EXCLUSION COLUMNS



This LC column is purpose-built for gene therapy and CAR-T workflows which rely heavily on lentiviral vectors (LVVs). Featuring an ultra-wide 2000 Å pore architecture engineered for the toughest large biomolecule separations, from the analysis of lentivirus aggregation to lipid nanoparticles.

TARGET APPLICATIONS


 Lentiviral Vectors Aggregate profiling & intact particle recovery	 Lipid Nanoparticles mRNA-LNP formulation & stability	 Plasmid DNA High MW nucleic acid separation	 sgRNA & Long Oligos CRISPR guide RNA characterization	 Viral-Like Particles VLP size distribution & purity	 PEG/Lipid Conjugates "Sticky" biologic analysis
---	--	---	---	---	---

KEY ADVANTAGES

2000 Å Ultra-Wide Pore Architecture Purpose-built for viral-sized particles. Achieves true size-exclusion of lentiviral vectors, aggregates, and LMW impurities without pore occlusion.	Hydrophilic Coated Surface Minimizes nonspecific interaction, resulting in improved recovery, enhanced peak symmetry, and more consistent size-based separation.
300 X Exceptional Column Lifetime Negligible performance deterioration after 300+ injections (150 mM phosphate, pH 7.0). Full-coverage bonded silica packing for unrivaled stability.	Versatility Detector Flexibility Low-noise baseline ensures seamless compatibility with UV, FLD, and SEC-MALS detection. Enables absolute determination of particle size and aggregation state.

TECHNICAL SPECIFICATIONS

Particle (µm)	5
Pore (Å)	2000
ID (mm)	4.6, 7.8
Length (mm)	50, 300
pH range	5.0-7.5 (optimal) 2.5-8.5 (extended)
Max pressure (psi)	2,000
Max temp (°C)	80
MP compatibility	Phosphate, TRIS, NH ₄ OAc, aqueous
Detection	UV, FLD, SEC-MALS



METHOD OPTIMIZATION INSIGHT
Ionic strength tuning is critical for lentivirus SEC. Moderate salt increases (e.g., 50 mM KCl) reduce electrostatic interactions for better HMW aggregate resolution, while excessive ionic strength reduces LMW selectivity. Optimal starting condition: 80 mM phosphate, 50 mM KCl, pH 6.5.

SURFACE TECHNOLOGY
Proprietary nanometer-thick hydrophilic film chemically bonded to high-purity spherical silica. The specialized neutral coating minimizes secondary interactions with lipid-enveloped viral particles, PEG-conjugated biologics, and other inherently "sticky" large biomolecules, ensuring separation is governed strictly by size exclusion.

Lentiviral vectors LNPs sgRNA Plasmid DNA
Water-soluble polymers VLPs PEG conjugates



COLUMN SCREENING SERVICE

Not sure which column or mobile phase to use? Sepax offers a full column screening and method development service. Rapid, cost-effective, and tailored to your sample type.

Scan for Order Information

Sepax Technologies, Inc.
Sepax-Tech.com
1-877-SEPAX-US
Info@Sepax-Tech.com

